



hekaDrill




ULTIMATE PRECISION

The hekaDrill has been developed with innovation and precision as a guiding principle.

Key capabilities include:

- High speed motor control - up to 85,000 rpm. Allows high torque/low speed control. (min. 200 rpm)
- Advanced motor control ensures comfort for the surgeon during extended use with minimal heat gain.
- Adjustable bur exposure lengths.
- Wireless foot pedal operates using a secure encrypted network, ensuring greater user protection.
- Adjustable acceleration and braking functionality.




A single rotary bur allows 5 possible exposure positions offering total flexibility and reduced inventory costs.

ELECTRIC HANDPIECE AND ATTACHMENTS

HIGHSPEED HANDPIECE		Order code	Model	Description	Speed rpm
		103347	HDS-121	Cable length 4m	200 - 85,000
		103352	HDS-122	Lever Control	200 - 85,000
		104184	Lever		
ATTACHMENTS STRAIGHT		Order code	Model	Length - mm	Bur size - mm
		103360	S-C07	70	2.38
		103371	S-C08	80	2.38
		103383	S-C10	100	2.38
		103394	S-C15	150	2.38
		103409	S-D09	90	3.18
		103414	S-D14	140	3.18
		103438	S-D26	260	3.18
SLIMLINE					
		104035	S-C07S	70	2.38
		104043	S-C08S	80	2.38
		104056	S-C10S	100	2.38
		104312	S-C12S	120	2.38
		104062	S-C15S	150	2.38
ANGLED					
		103446	A-C07	70	2.38
		103451	A-C10	100	2.38
		103467	A-C15	150	2.38
		103472	A-D14	140	3.18
SLIMLINE					
		104070	A-C07S	70	2.38
		104081	A-C10S	100	2.38
		104320	A-C12S	120	2.38
		104097	A-C15S	150	2.38
ANGLED MODULAR					
		104159	A-M01	Angled modular attachment	
CRANIOTOME					
		103508	C-C09	90 (Small)	2.38
		103480	C-C10	100 (Medium)	2.38
		103495	C-D12	120 (Large)	3.18
METAL CUTTER		Order code	Model	Speed rpm	Chuck type
		104137	M-D06	-	3.18
PERFORATOR ATTACHMENT					
		103519	P-P01	1,000	Hudson
		Used with HDS-121 at 72,000 rpm			
		103686	Locking spanner for HDS-1xx supplied with handpiece		
		Ask for catalogue 2677 for our range of rotary burs			



Description	Model	Order code
Console - less irrigation	PC-120	103305
Console - with irrigation	PC-121	103318



The advanced hekaDrill power console is available with or without irrigation and offers numerous features:

- Simultaneous dual handpiece operation.
- Adjustable acceleration and braking functionality.
- Integrated and customisable irrigation control with anti drip mode.
- 200+ user profiles can be stored and easily accessed.
- Rings illuminate to indicate port status for ease of use.

IRRIGATION AND FOOT PEDALS

STERILE IRRIGATION KIT



Sterile irrigation kit.
Box of 5.
Order code: S101461

FOOT PEDALS



- Both wired and wireless foot pedal available.
- Wireless foot pedal operates using a secure encrypted network for complete protection.
- A or B tool selection.
- Foot pedal underlighting to allow for easy detection.



4 metre cable for connecting and charging foot pedals via power console.

Order code: 103597

	Model	Order code
Wired	MFP-120	103556
Wireless	MFP-121	103562

POWER SUPPLY

Mains power cable 3 metres long. Please specify from below.

UK	Order code: 8290
Euro	Order code: 8300
USA/Canada	Order code: 8310
Australia/New Zealand	Order code: 8320
Japan	Order code: 8330
Swiss	Order code: 14700



HANDPIECE	Order code	Model	Speed rpm
-----------	------------	-------	-----------



103339	HDS-111	100,000
--------	---------	---------

The HDS-111 handpiece accepts all straight and angled attachments and craniotome

PERFORATOR	Order code	Model	Speed rpm	Operating pressure
------------	------------	-------	-----------	--------------------



103524	P-P02	850	7 bar
--------	-------	-----	-------

Uses Hudson chuck type

FOOT CONTROL	Order code	Model	Foot pedal inlet connection	Foot pedal outlet connection
--------------	------------	-------	-----------------------------	------------------------------



14900	FV-220	Integral 3 metre hose with MA7 probe	Accepts 'wash assist' hose from sizes below
-------	--------	--------------------------------------	---



3 metre length - Order code 632333
 4 metre length - Order code 632343
 5 metre length - Order code 632353



Wash cap for use with HDS-111/ P-P02 Order code 9950

STERILISATION TRAYS AND CONTAINER

STERILISATION TRAYS AND CONTAINER



ELECTRIC INSTRUMENT SYSTEM

Order code: 103610
Includes lid, container, electric tray and attachment tray.

Handpiece tray
Optional tray contents -
• 1 off HDS-121 handpiece
• 1 off P-P01 perforator
• 1 off electric cable
• 1 off spanner

Attachment tray
Optional tray contents -
• 5 off attachments
• 1 off spanner

PNEUMATIC INSTRUMENT SYSTEM

Order code 103629
Includes lid, container, pneumatic tray and attachment tray.

Handpiece tray
Optional tray contents -
• 1 off HDS-111 handpiece
• 1 off P-P02 perforator
• 1 off pneumatic hose
• 1 off spanner

Attachment tray
Optional tray contents -
• 5 off attachments
• 1 off spanner

CONTAINER AND TRAYS ONLY

Wire container and lid
Order code: 103944

Tray for electric system
Order code: 103953

Tray for pneumatic system
Order code: 103989

Tray for attachments
Order code: 103966

UNITED KINGDOM (INTERNATIONAL HEAD OFFICE)

De Soutter Medical Limited
Halton Brook Business Park, Weston Road,
Aston Clinton, Aylesbury, Bucks HP22 5WF
Tel: +44 (0) 1296 634 000
Email: info@de-soutter.com

AUSTRALIA

De Soutter Medical Australia Pty Ltd
2/12-14 Apollo Drive, Hallam
Victoria 3803
Tel: +61 (0) 3 9702 4441
Email: australia@de-soutter.com

AUSTRIA

De Soutter Medical Austria
Zweigniederlassung Österreich
Dietrichsteingasse 10, A-3400 Klosterneuburg
Tel: +43 (0) 676 96 71 770
Email: austria@de-soutter.com

BELGIUM

De Soutter Medical Belgium
Bessemmerstraat 14
3620 Lanaken
Tel: +32 (0) 89/47 15 37
Email: belgium@de-soutter.com

FRANCE

De Soutter Medical France
1 252 Avenue du Parc des Expositions
33260 La Teste de Buch
Tel: +33 (0) 5 56 54 89 36
Email: france@de-soutter.com

GERMANY

De Soutter Medical Deutschland
Bahnhofstraße 4
D-66625 Nohfelden
Tel: +49 (0) 68 52-99 12 46
Email: deutschland@de-soutter.com

ITALY

De Soutter Medical Italia
Località Fornace SNC
27022 Casorate Primo - PV
Tel: +39 (0) 2 9009 4098
Email: italy@de-soutter.com

NETHERLANDS

De Soutter Medical Nederland
Gelderlandhaven 2X
3433 PG Nieuwegein
Tel: +31 (0) 85 0491480
Email: nederland@de-soutter.com

UNITED STATES OF AMERICA

De Soutter Medical USA, Inc.
224 Rolling Hill Road, Suite 12A
 Mooresville, NC 28117
Tel: +1 (704) 655 9040
Email: usa@de-soutter.com

deSoutter
MEDICAL

Distributed by De Soutter Medical, manufactured by Zethon